KORISNIČKI PRIRUČNIK

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- Detaljno objašnjenje svake ikone nalazi se na posljednjoj stranici.

LG DXD Calibration

Calibration User BPM Validation FI

DXD Calibration SW Firmware Update About

Connection & File Save

DXD Serial Number

DXD IP . . .

Timeout msec, (3000-10000)

Save Location

Network Options

Current Status: Wired

Installation Info.

Date Format

Current Date 2018/01/19

Detector Parameters

Trigger Mode ☐ Auto ☒ Manual

Sensitivity (0-63)

Window Time 00 msec (1-40)

Frame Width

Frame Height

! Press 'Reset' to load factory-default Detector

Power Options

Auto Sleep

Auto Power-Off

Date	Time	Type	Details
2018-01-19	15:04:12	Settings	Ping succeeded to 10.10.10.100
2018-01-19	15:04:14	Settings	Failed to load factory calibration files
2018-01-19	15:04:14	Settings	Connected to detector successfully

- [Calibration SW]: konfiguracija postavki potrebnih za kalibraciju algoritama Calibration Software.

! NAPOMENA

- Detaljno objašnjenje svake ikone nalazi se na posljednjoj stranici.

LG DXD Calibration

Calibration
User BPM
Validation
E

DXD
Calibration SW
Firmware Update
About

Calibration Parameters

Target Gain
(0.0~255.0)

Gain Margin
(0.0~1.0)

Offset Margin
(0~1000)

Std Margin
(0~1000)

Ref sat value
(28000~65534)

Surr Margin
(0~3000)

Cut Edge

T
(6~1000)

L
(6~1000)

R
(6~1000)

B
(6~1000)

Image Edit

Rotation

Flip

Invert
☐

Press 'Reset' to load factory-default Calibration SW

Reset

Save
Cancel

Exit

Date	Time	Type	Details
2018-10-01	09:22:18	Settings	Ping succeeded to 10.10.10.100
2018-10-01	09:22:21	Settings	Connected to detector successfully

- [Firmware Update]: provjera verzije firmvera detektora ili ažuriranje firmvera. U ovom izborniku možete ažurirati firmver.

The screenshot shows the 'LG DXD Calibration' application window. The 'Firmware Update' tab is selected, displaying the 'Current Firmware' section with 'Version 3.00.01'. The 'New Firmware' section includes a 'Firmware File' input field with an 'Open' button, a 'Status' label, and 'Update' and 'Factory Reset' buttons. An 'Exit' button is located at the bottom right. A log table at the bottom shows two entries from 2019-09-26.

Date	Time	Type	Details
2019-09-26	13:53:42	Settings	Ping succeeded to 10.10.10.100
2019-09-26	13:53:45	Settings	Connected to detector successfully

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

The screenshot displays the 'LG DXD Calibration' software window with the 'Validation' tab selected. The interface includes a 'Navigator' on the left, a central 'Pixel Value' display, a 'Histogram' with a 'Reset' button, and a 'History' panel on the right. The 'BPM Analysis' table is visible in the bottom left, and a log of events is shown at the bottom.

BPM Analysis

Item	Count
Rows	0
Columns	0
Class 0	0
Class 1	0
Class 2	0
Class 3~4	0, 0
Class 5~6	0, 0
Class 7~8	0, 0

History

Dark Average Map ☐
Offset Map ☐
Gain Map ☐
Bad Pixel Map ☐

Image(s) Count: 0 Acquisition

Image not Loaded

Date	Time	Type	Details
2018-10-01	09:22:18	Settings	Ping succeeded to 10.10.10.100
2018-10-01	09:22:21	Settings	Connected to detector successfully

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
History

Image(s)	Count: 0	Acquisition
	Dose	Median
Measure & Save		

Image not Loaded

Date	Time	Type	Details
2018-10-01	09:22:18	Settings	Ping succeeded to 10.10.10.100
2018-10-01	09:22:21	Settings	Connected to detector successfully

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Result

Image not Loaded

Date	Time	Type	Details
2018-10-01	09:22:18	Settings	Ping succeeded to 10.10.10.100
2018-10-01	09:22:21	Settings	Connected to detector successfully

Pokretanje programa

- Dvaput kliknite na izvršnu datoteku koju ste instalirali na računalo za pokretanje Calibration Software.
- Prilikom prvog pokretanja prikazat će se zaslon s postavkama.

The screenshot shows the 'LG DXD Calibration' application window. It has a menu bar with 'Calibration', 'User BPM', 'Validation', and 'EI'. Below the menu bar are four tabs: 'DXD' (selected), 'Calibration SW', 'Firmware Update', and 'About'. The 'DXD' tab contains several sections:


- Connection & File Save:** Includes fields for 'DXD Serial Number', 'DXD IP' (10.10.10.100), 'Timeout' (5000 msec), and 'Save Location' (D:\). There are 'Ping' and 'Open' buttons.
- Network Options:** Shows 'Current Status' as 'Wired' with buttons for 'Change DXD IP' and 'Wireless Setup'.
- Installation Info:** Includes 'Date Format' (YYYY/MM/DD) and 'Current Date' (2018/01/19) with a 'Register' button.
- Detector Parameters:** Includes 'Trigger Mode' (Auto/Manual), 'Sensitivity' (15), 'Window Time' (5.00 msec), 'Frame Width' (2500), and 'Frame Height' (3052). A warning message says 'Press "Reset" to load factory-default Detector' with a 'Reset' button.
- Power Options:** Includes 'Auto Sleep' and 'Auto Power-Off' (both set to Off) with 'Save' and 'Cancel' buttons.

At the bottom of the window is a table with the following data:


Date	Time	Type	Details
2018-01-19	15:04:12	Settings	Ping succeeded to 10.10.10.100
2018-01-19	15:04:14	Settings	Failed to load factory calibration files
2018-01-19	15:04:14	Settings	Connected to detector successfully

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The screenshot shows the 'LG DXD Calibration' application window. It has a dark-themed header with tabs: Calibration, User BPM, Validation, and EI. Below the header, there are sub-tabs: DXD, Calibration SW, Firmware Update, and About. The 'DXD' sub-tab is active, displaying several configuration sections:

- Connection & File Save:** Includes fields for DXD Serial Number, DXD IP (10.10.10.100), Timeout (5000 msec), and Save Location (D:\). There are 'Ping' and 'Open' buttons.
- Network Options:** Shows 'Current Status' as 'Wired' with buttons for 'Change DXD IP' and 'Wireless Setup'.
- Installation Info:** Includes 'Date Format' (YYYY/MM/DD) and 'Current Date' (2018/01/19) with a 'Register' button.
- Detector Parameters:** Includes 'Trigger Mode' (Auto/Manual), 'Sensitivity' (15), 'Window Time' (5 00 msec), 'Frame Width' (2500), and 'Frame Height' (3052). A warning message says 'Press 'Reset' to load factory-default Detector'.
- Power Options:** Includes 'Auto Sleep' and 'Auto Power-Off', both set to 'Off'.

At the bottom, there is a log table with columns: Date, Time, Type, and Details.

Date	Time	Type	Details
2018-01-19	15:04:12	Settings	Ping succeeded to 10.10.10.100
2018-01-19	15:04:14	Settings	Failed to load factory calibration files
2018-01-19	15:04:14	Settings	Connected to detector successfully

- [Detector Parameters]: postavke koje se upotrebljavaju prilikom dobivanja snimke putem detektora.
- Kliknite na gumb [Save] kako biste primijenili unesene postavke.
- Detalji postavki navedeni su u nastavku:
 - [Trigger Mode]: postavljanje načina aktivacije [Trigger Mode].
 - * [Auto]: omogućavanje značajke automatskog otkrivanja izlaganja.
 - * [Manual]: onemogućavanje značajke automatskog otkrivanja izlaganja.
 - [Sensitivity]: osjetljivost ploče.
 - [Window Time]: postavljanje vremena za čitanje podataka nakon izlaganja rendgenskom zračenju (jedinica: 100 ms; kada se unese 5, vrijeme se postavlja na 500 ms).
 - [Frame Width] / [Frame Height]: broj piksela detektora.
- Radnja svakog gumba opisana je u nastavku:
 - [Save]: primjena podešenih postavki.
 - [Reset]: učitavanje tvorničkih postavki.
 - [Cancel]: učitavanje posljednjih spremljenih postavki.

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The screenshot shows the 'LG DXD Calibration' software window. The 'Calibration SW' tab is active, displaying 'Calibration Parameters' and 'Image Edit' sections. The 'Calibration Parameters' section includes input fields for Target Gain (1), Gain Margin (0.6), Offset Margin (120), Std Margin (80), Ref sat value (56000), Surr Margin (2000), and Cut Edge (T: 6, L: 6, R: 6, B: 6). The 'Image Edit' section includes dropdowns for Rotation (0) and Flip (None), and a checkbox for Invert. Below these sections are buttons for 'Reset', 'Save', 'Cancel', and 'Exit'. At the bottom, a log table shows the following data:

Date	Time	Type	Details
2018-10-01	09:22:18	Settings	Ping succeeded to 10.10.10.100
2018-10-01	09:22:21	Settings	Connected to detector successfully

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


Date	Time	Type	Details
2018-10-01	09:53:49	User BPM	Image saved
2018-10-01	09:53:50	User BPM	Image opened
2018-10-01	09:53:50	User BPM	Acquisition end

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

Date	Time	Type	Details
2018-10-01	10:15:53	User BPM	[FILE INFO] avg_dark_size=18874368, offset_size=18874368, gain_size=37748736, bpm_size=18874368
2018-10-01	10:15:53	User BPM	[Warning] Do not click any buttons or tab during upload.
2018-10-01	10:15:56	User BPM	Package file creation success
2018-10-01	10:15:56	User BPM	Uploading calibration map to dxd... map size = 94371872
2018-10-01	10:15:56	User BPM	[Warning] Do not click any buttons or tab during upload.
2018-10-01	10:16:02	User BPM	New map upload success. Previous map files are all deleted.

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Gain Map ☒

Bad Pixel Map ☒


Date	Time	Type	Details
2018-10-01	10:11:01	Validation	Image saved
2018-10-01	10:11:02	Validation	Image opened
2018-10-01	10:11:02	Validation	Acquisition end

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LG DXD Calibration

Calibration User BPM Validation EI

DXD Calibration SW Firmware Update **About**

Connection & File Save

DXD Serial Number 12345678

DXD IP 10 . 10 . 10 . 100

Timeout 5000 msec. (3000~10000) Ping

Save Location D:\ Open

Apply

Network Options

Current Status: Wired Change DXD IP Wireless Setup

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




Apply

Network Options

Current Status: Wired **Change DXD IP** Wireless Setup

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

[Apply]

Network Options

Current Status: Wired [Change DXD IP] [**Wireless Setup**]

Installation Info.

Date Format: [YYYY/MM/DD]

Current Date: 2018/12/13 [Register]

Detector Parameters

Trigger Mode: ☒ Auto ☐ Manual

Sensitivity: [15] (0-63)

Window Time: [5] 00 msec (1-40)

Frame Width: [2500]

Frame Height: [3052]

! Press 'Reset' to load factory-default Detector [Reset]

[Save] [Cancel]

Power Options

Auto Sleep: [Off]

Auto Power-Off: [Off]

[Save] [Cancel]

[Exit]

Date	Time	Type	Details
2018-12-13	13:33:48	Settings	Settings: Invalid factory calibration files
2018-12-13	13:33:49	Settings	Connected to detector successfully
2018-12-13	13:34:05	Settings	Restored factory-default parameters

- Ako se prikaže skočni prozor, unesite svoj [SSID] i [Password] i zatim kliknite [Apply].

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







DXD Monitoring System			
Product Information	Network Information	Battery	ETC
Software Version Firmware Ver. 2.00.01	Status  Wired Connected	Status  Battery Connected  Fully Charged 100 %	Bright Image Count 1
Manufacturing Date 2017. 11. 20	IP 10.10.10.100	Auto Sleep Off	Dark Image Count 4
Installation Date 0000. 00. 00	SSID N/A	Auto Power-Off Off	
Model Number 14HK701G	Netmask 255.255.255.0		
Serial Number 12345678	Gateway 10.10.10.1		
	Mac 78:5D:C8:B9:44:09		

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



- 1 Idite na [Control Panel] i odaberite izbornik [System and Security].

- 2 Kliknite na poveznicu [Windows Firewall].


The screenshot shows the Windows Action Center interface. It has three main sections: 'Action Center' with a flag icon, 'Windows Firewall' with a globe and brick wall icon, and 'System' with a monitor icon. The 'Windows Firewall' section is highlighted with a light blue background. Within this section, the link 'Allow a program through Windows Firewall' is highlighted in blue. Other links include 'Review your computer's status and resolve issues', 'Change User Account Control settings', 'Troubleshoot common computer problems', 'Restore your computer to an earlier time', 'Check firewall status', 'View amount of RAM and processor speed', 'Check the Windows Experience Index', 'Allow remote access', 'See the name of this computer', and 'Device Manager'.

Action Center
Review your computer's status and resolve issues | Change User Account Control settings |
Troubleshoot common computer problems | Restore your computer to an earlier time

Windows Firewall
Check firewall status | Allow a program through Windows Firewall

System
View amount of RAM and processor speed | Check the Windows Experience Index |
Allow remote access | See the name of this computer | Device Manager

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The screenshot shows a list of links in the Windows Firewall control panel. The link 'Advanced settings' is highlighted in blue. The list includes 'Control Panel Home', 'Allow a program or feature through Windows Firewall', 'Change notification settings', 'Turn Windows Firewall on or off', 'Restore defaults', 'Advanced settings', and 'Troubleshoot my network'.

Control Panel Home

Allow a program or feature through Windows Firewall

Change notification settings


Turn Windows Firewall on or off

Restore defaults

Advanced settings

Troubleshoot my network

- 4 Odaberite stavku [Inbound Rules] u opciji Windows Firewall with Advanced Security.

- 5 Pomaknite prema dolje i pronađite pravilo [File and Printer Sharing] te kliknite na [Enable Rule].

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
File and Printer Sharing (Echo Request - ICMPv4-In)	File and Printer Sharing	Private, Public
File and Printer Sharing (Echo Request - ICMPv4-In)	File and Printer Sharing	Domain
File and Printer Sharing (Echo Request - ICMPv6-In)	File and Printer Sharing	Domain
File and Printer Sharing (Echo Request - ICMPv6-In)	File and Printer Sharing	Private, Public

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
Wireless

View and change router settings

Wireless **MAC Filtering** Wi-Fi Protected Setup

Smart Connect: ☒ 5 GHz band steering ☐ Off
5 GHz + 5 GHz Individual networks

2.4 GHz Network: **ON**

Network name: Password: Broadcast SSID: Channel:

Security mode: Network mode: Channel width:

5 GHz + 5 GHz Network: **ON**

Network name: Password: Broadcast SSID: Channel:

Security mode: Network mode: Channel width:

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Copyright © All Rights Reserved