KORISNIČKI PRIRUČNIK

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Result

Image not Loaded

Date	Time	Type	Details
2018-10-01	09:22:18	Settings	Ping succeeded to 10.10.10.100
2018-10-01	09:22:21	Settings	Connected to detector successfully

RAD

Pokretanje programa

- Duplim klikom na datoteku koja se izvršava i koja je instalirana na računaru pokrenite Calibration Software.
- Kada je pokrećete prvi put, bićete preusmereni na ekran Podešavanja.

LG DXD Calibration

Calibration User BPM Validation EI

DXD Calibration SW Firmware Update About

Connection & File Save

DXD Serial Number

DXD IP . . .

Timeout msec (3000~10000)

Save Location

Network Options

Current Status: Wired

Installation Info.

Date Format

Current Date 2018/01/19

Detector Parameters

Trigger Mode ☐ Auto ☒ Manual

Sensitivity (0~63)

Window Time 00 msec (1~40)

Frame Width

Frame Height

! Press 'Reset' to load factory-default Detector

Power Options


Auto Sleep

Auto Power-Off


Date	Time	Type	Details
2018-01-19	15:04:12	Settings	Ping succeeded to 10.10.10.100
2018-01-19	15:04:14	Settings	Failed to load factory calibration files
2018-01-19	15:04:14	Settings	Connected to detector successfully

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The screenshot displays the 'LG DXD Calibration' application window. It features a top navigation bar with tabs: Calibration, User BPM, Validation, and EI. Below this is a sub-navigation bar with 'DXD', 'Calibration SW', 'Firmware Update', and 'About'. The main content area is divided into several sections:

- Connection & File Save:** Includes fields for DXD Serial Number, DXD IP (10.10.10.100), Timeout (5000 msec), and Save Location (D:\). Buttons for 'Ping' and 'Open' are present.
- Network Options:** Shows 'Current Status' as 'Wired' with buttons for 'Change DXD IP' and 'Wireless Setup'.
- Installation Info:** Includes 'Date Format' (YYYY/MM/DD) and 'Current Date' (2018/01/19) with a 'Register' button.
- Detector Parameters:** Includes 'Trigger Mode' (Auto/Manual), 'Sensitivity' (15), 'Window Time' (5 00 msec), 'Frame Width' (2500), and 'Frame Height' (3052). A warning message states: 'Press 'Reset' to load factory-default Detector'.
- Power Options:** Includes 'Auto Sleep' and 'Auto Power-Off', both set to 'Off'.

At the bottom, there is a log table with the following data:

Date	Time	Type	Details
2018-01-19	15:04:12	Settings	Ping succeeded to 10.10.10.100
2018-01-19	15:04:14	Settings	Failed to load factory calibration files
2018-01-19	15:04:14	Settings	Connected to detector successfully

- [Detector Parameters]: Podešavanja koja su korišćena prilikom preuzimanja snimka sa detektora.
- Kliknite na dugme [Save] da primenite uneta podešavanja.
- Detalji o podešavanjima su sledeći:
 - [Trigger Mode]: Podesite [Trigger Mode].
 - * [Auto]: Omogućite funkciju automatske detekcije izlaganja.
 - * [Manual]: Onemogućite funkciju automatske detekcije izlaganja.
 - [Sensitivity]: Osetljivost ploče.
 - [Window Time]: Podesite vreme za očitavanje podataka nakon izlaganja rendgenskim zracima.
(Jedinica: 100 ms, kada unesete 5, vreme se podešava na 500 ms)
 - [Frame Width] / [Frame Height]: Broj piksela na detektoru.
- Funkcije svakog dugmeta su sledeće:
 - [Save]: Primena promenjenih podešavanja.
 - [Reset]: Učitavanje fabričkih podešavanja.
 - [Cancel]: Učitavanje poslednjih očitanih podešavanja.

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





The screenshot shows the LG DXD Calibration software interface. The 'User BPM' tab is selected. The main image area displays a dark image with a yellow circle highlighting a pixel. The histogram shows a single peak. The BPM analysis table shows the following data:

Item	Count
Rows	0
Columns	0
Class 0	119
Class 1	18
Class 2	0
Class 3~4	8, 0
Class 5~6	0, 0
Class 7~8	0, 0

The history panel shows the following data:

Image(s)	Count	Acquisition
BPM_20181001095349.raw	1	

The pixel view table shows the following data:

1502	1503	1504	1505	1506	1507	1508
1520	12	10	9	13	14	9
1521	12	13	14	15	14	13
1522	12	21	18	12	14	16
1523	11	14	14	12	11	13
1524	13	10	15	10	14	11
1525	9	8	16	15	8	12


The bottom status bar shows the following information: Image Loaded, Min=0, Max=5434, Med=8, Mean=10.87, Std=82.07, Bad Pixel Count=145, (2488,2688) 16B=7, 8B=0, 12.50%, 3072x3072.

1 Uđite u meni [User BPM]

- Kliknite na meni [User BPM] da uđete.
- [User BPM] zahteva preuzimanje snimka jer vizuelno pregleda snimak na koji se rezultat [Calibration] primenjuje.

! NAPOMENA

- Sledeći iskačuci prozor se pojavljuje prilikom ulaska u meni bez obavljanja [Calibration].

2 Proverite datoteku sa istorijom


- Proverite da li se naziv [History] kreiran na osnovu [Calibration] poklapa sa nazivom prikazanim u prozoru aktuelne [History].
- Primenite izabranu datoteku Istorije i sprovedite postupak [User BPM].

3 Preuzimanje snimaka


- Kliknite na dugme [Acquisition] i preuzmite svetao snimak. Naziv snimka će se prikazati na listi [Image(s)].
- Informacije o snimku prikazuju se u prikazu snimka u nastavku.

4 [Pixel View]

- Proverite vrednosti piksela u [Pixel View].
 - Vrednosti piksela iz središta prikazivača snimka priloženi su u [Pixel View].
 - Tu su minimalna vrednost, maksimalna vrednost i mogući loši pikseli prikazani na sledeći način:
 - Minimalna vrednost: Prikazana plavim brojkama.
 - Maksimalna vrednost: Prikazana crvenim brojkama.
 - Mogući loši pikseli: Prikazani na sivoj pozadini.
 - Veličina prozora [Pixel View] može da se promeni korišćenjem ikonice . Ikonica  se pojavljuje kada se mišem prelazi preko granice između [Pixel View] i prikazivača snimaka.

- Podesite dodatne loše piksele u [Pixel View]
 - Dvostrukim klikom na piksel u [Pixel View] podesite piksel kao loš piksel. Ponovnim dvostrukim klikom na isti piksel koji je naznačen kao loš piksel otkažite izbor.
 - Ako je piksel podešen kao loš piksel, vrednost se ažurira u [BPM Analysis]. Naznačeni loš piksel se zamenjuje kalibrisanom vrednošću piksela.
- Sačuvajte konačni [User BPM].
 - Prilikom ulaska u drugi meni sačuvaće se datoteka sa rezultatom.
 - Iskašuća poruka se pojavljuje i od vas se traži da izaberete da li želite da sačuvate datoteku prilikom izlaska iz menija u kom ste trenutno i ulaska u drugi meni.
 - Kada se sačuva, jedna [History] se dodaje, a datoteka [BPM.raw] se ažurira i čuva.

Dodeljivanje loših linija u prikazu piksela

- Kada se naznačava loš piksel, to je funkcija za naznačavanje jedinice linije umesto jedinice piksela.
- Nakon što potvrdite polja za potvrđivanje za vertikalni ili horizontalni smer, kliknite dva puta na piksel u [Pixel View] da naznačite liniju u izabranom smeru.
 - Ako naznačite liniju od tačke 1.522 do tačke 2.600, izaberite dugme [Apply] nakon što unesete vrednost da biste dodelili liniju kao lošu liniju.

Pixel View ☐ H-Bad Line ☒ V-Bad Line

	1502	1503	1504	1505	1506	1507	1508
1520	12	10	9	13	14	13	9
1521	12	13	14	15	14	13	14
1522	12	21	18	14	14	16	16
1523	11	14	14	12	11	13	10
1524	13	10	15	10	14	12	11
1525	9	8	16	15	8	12	12

Pixel (x 1505, y 1522) has been updated.
Bad Pixel Map Upload

Pixel View ☐ H-Bad Line ☐ V-Bad Line

	1502	1503	1504	1505	1506	1507	1508
1520	12	10	9	13	14	13	9
1521	12	13	14	15	14	13	14
1522	12	21	18	14	14	16	16
1523	11	14	14	12	11	13	10
1524	13	10	15	10	14	12	11
1525	9	8	16	15	8	12	12

Pixel (x 1505, y 1522) has been updated.
Bad Pixel Map Upload

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Date	Time	Type	Details
2018-10-01	10:15:53	User BPM	[FILE INFO] avg_dark_size=18874368, offset_size=18874368, gain_size=37748736, bpm_size=18874368
2018-10-01	10:15:53	User BPM	[Warning] Do not click any buttons or tab during upload.
2018-10-01	10:15:56	User BPM	Package file creation success
2018-10-01	10:15:56	User BPM	Uploading calibration map to dxd... map size = 94371872
2018-10-01	10:15:56	User BPM	[Warning] Do not click any buttons or tab during upload.
2018-10-01	10:16:02	User BPM	New map upload success. Previous map files are all deleted.

! NAPOMENA

- [Bad Pixel Map] Prilikom otpremanja, izbrisite mapu u postojećem folderu Fabrička kalibracija. Da biste je sačuvali, to se mora učiniti nakon pravljenja rezervne kopije na zasebnoj putanji.
- [Bad Pixel Map] Prilikom biranja datoteke za otpremanje, sve četiri datoteke [Bad Pixel Map] moraju da budu na putanji. ([BPM.raw], [AvgDark.raw], [Offset.raw], [Gain.raw])

Validation

Ovaj meni omogućava korisnicima da vizuelno provere rezultat [Calibration] nakon obavljanja [Calibration].

1

2

3

4

LG DXD Calibration

Calibration User BPM Validation E

Navigator + - 1530 1530 Pixel Value: 15

Histogram Reset

Window 65536

Level 32768

BPM Analysis

Item	Count
Rows	0
Columns	0
Class 0	119
Class 1	18
Class 2	0
Class 3~4	8, 0
Class 5~6	0, 0
Class 7~8	0, 0

Image Loaded Min=0, Max=4594, Med=7, Mean=10.78, Std=87.62, Bad Pixel Count=145 100.00% 3060x3060

History Open

Factory Calibration

Images(s) Count: 1 Acquisition

Validation_20181001101101.raw

Dark Average Map ☒

Offset Map ☒

Gain Map ☒

Bad Pixel Map ☒


Date	Time	Type	Details
2018-10-01	10:11:01	Validation	Image saved
2018-10-01	10:11:02	Validation	Image opened
2018-10-01	10:11:02	Validation	Acquisition end

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





The screenshot shows the 'LG DXD Calibration' application window. The 'About' button is located in the top right corner of the main window, highlighted with a red rectangle. The application has a dark theme and a sidebar on the left with tabs: Calibration, User BPM, Validation, and EI. The main area is divided into several sections: Connection & File Save, Network Options, Installation Info., Detector Parameters, and Power Options. The 'About' button is located in the top right corner of the main window, highlighted with a red rectangle.

Connection & File Save

DXD Serial Number: 12345678

DXD IP: 10 . 10 . 10 . 100

Timeout: 5000 msec. (3000~10000) Ping

Save Location: D:\ Open

Network Options

Current Status: Wired Change DXD IP Wireless Setup

Installation Info.

Date Format: YYYY/MM/DD

Current Date: 2018/10/01 Register

Detector Parameters

Trigger Mode: ☐ Auto ☒ Manual

Sensitivity: 15 (0~63)

Window Time: 5 00 msec (1~40)

Frame Width: 3072

Frame Height: 3072

! Press 'Reset' to load factory-default Detector Reset

Power Options

Auto Sleep: Off

Auto Power-Off: Off

Save Cancel

Exit

Date	Time	Type	Details
2018-10-01	09:22:18	Settings	Ping succeeded to 10.10.10.100
2018-10-01	09:22:21	Settings	Connected to detector successfully

The 'About' dialog box is shown. It has a title bar with the text 'About' and a close button. The main content area contains the following text:

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






The screenshot shows the 'LG DXD Calibration' software window. The 'Calibration' tab is active, and the 'Wireless Setup' button is highlighted with a red rectangle. The interface is divided into several sections:

- Connection & File Save:** Includes fields for DXD Serial Number, DXD IP (10.10.10.100), Timeout (5000 msec), and Save Location (D:\).
- Detector Parameters:** Includes Trigger Mode (Auto), Sensitivity (15), Window Time (500 msec), Frame Width (2500), and Frame Height (3052).
- Network Options:** Includes Current Status (Wired) and a 'Wireless Setup' button.
- Installation Info:** Includes Date Format (YYYY/MM/DD) and Current Date (2018/12/13).
- Power Options:** Includes Auto Sleep (Off) and Auto Power-Off (Off).

A log table at the bottom shows the following entries:

Date	Time	Type	Details
2018-12-13	13:33:48	Settings	Settings: Invalid factory calibration files
2018-12-13	13:33:49	Settings	Connected to detector successfully
2018-12-13	13:34:05	Settings	Restored factory-default parameters

- Ako se pojavi iskačući prozor, unesite [SSID] i [Password], i kliknite na [Apply].

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Wireless Setup [X]

Wireless Setup is complete.
Perform Ping Test again.

*If you connect to the new IP address,
the existing connection will be disconnected.

<Iskačući prozor za uspešnu konfiguraciju>


Wireless Setup [X]

Cannot complete Wireless Setup.
Please try again.

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






DXD Monitoring System			
Product Information	Network Information	Battery	ETC
Software Version Firmware Ver. 1.00.18	Status  Wired Connected	Status  Battery Connected  Fully Charged 100 % 	Bright Image Count 52
Manufacturing Date 2017. 06. 03	IP 10.10.10.100	Auto Sleep Off	Dark Image Count 28
Installation Date 0000. 00. 00	SSID N/A	Auto Power-Off Off	
Model Number 17HK701G	Netmask 255.255.255.0		
Serial Number 12345678	Gateway 10.10.10.1		
	Mac 0E:74:B1:23:EA:85		

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


- 2 Kliknite na vezu [Windows Firewall].

Action Center

[Review your computer's status and resolve issues](#) | [Change User Account Control settings](#) |
[Troubleshoot common computer problems](#) | [Restore your computer to an earlier time](#)

Windows Firewall

[Check firewall status](#) | [Allow a program through Windows Firewall](#)

System

[View amount of RAM and processor speed](#) | [Check the Windows Experience Index](#) |
[Allow remote access](#) | [See the name of this computer](#) | [Device Manager](#)

- 3 Na levoj strani prozora kliknite na vezu [Advanced Settings].

Control Panel Home

[Allow a program or feature through Windows Firewall](#)

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



Wireless

View and change router settings

Wireless MAC Filtering Wi-Fi Protected Setup

Smart Connect: ☒ 5 GHz band steering ☐ Off
5 GHz + 5 GHz Individual networks

Network name: **LGEDXD** 2.4 GHz Network: ☒ ON

Password: **[REDACTED]** Broadcast SSID: Yes Channel: Auto

Security mode: WPA2 Personal Network mode: Mixed Channel width: Auto

Network name: **LGEDXD** 5 GHz + 5 GHz Network: ☒ ON

Password: **[REDACTED]** Broadcast SSID: Yes Channel: Auto

Security mode: WPA2 Personal Network mode: Auto Channel width: Auto

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Copyright © All Rights Reserved