LG Pra.L

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Quick Guide

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Facial Contouring Device

BLQ1

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ENGLISH

Please read this manual carefully before operating your set and retain it for future reference.

Product Features


Improving skin elasticity with the micro-focus ultrasound

The product focuses ultrasound energy on the target area of a skin layer, creating a point of energy coagulation and promoting the generation of new collagen and elastin in the area.


※ The produced effects may vary depending on the skin properties of the user and the cosmetics that the user uses.

※ This product is a skin care device, not a medical device, and produces no medical effect whatsoever other than the cosmetic effects related to the product function.

Components


Exploring the Product

- 1 Ultrasound irradiation part
- 2 Cartridge mount/unmount button
- 3 Status indicator
- 4 LEVEL (intensity) button
- 5 SHOT button
- 6 Power button
- 7 Strap hole
- 8 Speaker
- 9 Charging port (C type)
- 10 Vent
- 11 Charging contact

Charging

Product Specifications

General	
Dimensions (W x H x D)	Main body: Approx. 36 mm x 40 mm x 200 mm Charging cradle: Approx. 50.8 mm x 56.2 mm x 235.7 mm
Operating temperature	5°C - 35°C
Operating humidity	0% - 90%
Storage temperature	1 °C - 60 °C ※ If the storage temperature is 0 °C or lower, the cartridge may freeze and burst, so be careful.
Storage humidity	0% - 90%
Recommended number of uses	Once a week, 300 shots each time • When fully charged, you can use the product four times (300 shots each) at LEVEL 5. (Internal test result) • The hours of use and lifetime can vary depending on the battery warranty period, ambient temperature and conditions of use.
Time to full charge	Less than 4 hours • Based on a 5 V, 2.0 A charging adapter • Charging time may vary depending on charging adapter specifications, battery level, ambient temperature and charging conditions.

• Design and specifications are subject to change without notice.

Precautions for Use

Before performing care

- Those who have a facial skin disease or frequently suffer from facial dermatitis should consult a dermatologist before use.
- As a product used on the face and neck, it is not recommended to use on other parts of the body.
- Do not use the product over the areas injected with Botox, fillers or autologous fat or the areas treated with thread lift.
- If crack or damage occurs on the head part that touches the skin, stop using the product and receive service at the service center.

During care

- Apply moisture gel on the device's irradiation part or guide part thoroughly before use.
- When using the product, put the applicator against the skin tightly.
- Do not apply excessive force when contacting the product to skin. When pushing with excessive force, a point of energy coagulation is created in an area deeper than the skin layer, which may cause more severe pain or other troubles. Therefore, adhere to the skin with appropriate force when using the product.
- Do not irradiate in the air without appropriate skin contact.
- Avoid care between the eyebrows, around the eyes, around the mouth (1cm from the oral angle), jawbone line, and uvula line.
- Adjust the intensity according to the skin properties and do not fire shots to the same area repeatedly or continuously.
- When you use the product for the first time, start with LEVEL (intensity) 1. Increase the level as the number of uses increases.
- Using the product more frequently than the recommended use frequency may cause skin irritation. (Once a week, total 300 shots in the guide parts)

After care

- When this product is used, there may be skin irritation such as blushing or tingling depending on the user's skin condition. If such irritation persists, stop using it immediately and consult a dermatologist before using it again.
- Avoid places where skin may be overstimulated, such as sauna or spa, on the day of treatment.
- Avoid using scrub or peeling agent for about 2-3 days after the treatment.

Safety cautions

- CAUTION concerning the Power Cord**
The Power Plug is the disconnecting device. In case of an emergency, the Power Plug must remain readily accessible.
- CAUTION:** The DC Adapter must be purchased separately, and a certified DC adapter (5 V, 2.0 A or higher) must be used. If you use an uncertified charger, there is a risk of breakdown and change in charging time. If charging is unsuccessful, it is recommended to try again with an alternative certified DC adapter.
- When charging the product, be cautious not to get the charging part to be stained with foreign substances such as moisture or cosmetics. Otherwise, it may cause damage to the product or a fire.

For safety precautions, refer to the manual provided separately.

Q&A

Q How do I maintain the product?


Do not wash with water as it is not waterproof. It may cause product malfunction. If there are foreign substances left on the product or product head after the care, wipe off that part with water tissue and wipe again with a soft dry towel or cloth.

Do not use strong chemicals such as alcohol, benzene, or acetone.

Chemical products (alcohol, benzene, acetone, etc.)

- Make sure that no water or cosmetics gets inside the vent or gap in the product.
- Make sure that the film of the irradiation part is not torn by a sharp object such as a knife or nail.
- This product is vulnerable to shock, so make sure that it does not fall to the ground and gets damaged.

Getting to Know the Status Indicator


- 1 Number of shots fired / Total number of shots recommended
- 2 Voice guide on / Voice guide off (Mute)
- 3 Battery level / Battery charging status
- 4 User number

Product status display


- READY:** The product power is turned on and ready for use.
- Number:** Indicates the number of shots fired.
- REPLACE:** There is no shot left in the cartridge. Please replace the cartridge.
- NO CARTRIDGE:** The cartridge is absent or not mounted properly. Please mount the cartridge properly.
- NO BATTERY:** The battery has exhausted. Please recharge the battery.
- TILT:** The product is tilted and deviates from the recommended use angle. Please adhere to the recommended use angle. (Please refer to "How to Use the Device Correctly" in the Care Guide.)
- CHARGING...:** The product is turned on while being charged. Please use the product after charging is complete.
- Care LEVEL (intensity): 1 - 5

Setting Voice Guide


- Voice guide on → Sound effect → Voice guide off (Mute)
- When using for the first time, follow the voice guide. Once you get familiarized with appropriate use, change to the sound effect or mute and use the product.

Changing the User


While the power of the product is turned on by pressing and holding the Power button for 2 seconds or longer, every time you press the LEVEL (intensity) button on the front and the Power button on the back for about 3 seconds at the same time, a voice guide setting changes.

Setting Language

While the power of the product is turned on by pressing and holding the Power button for 2 seconds or longer, every time you press the LEVEL (intensity) button on the front and the SHOT button on the back for about 3 seconds at the same time, a used language changes.

Removing/Reconnecting the Strap from the Main Body

The strap is provided connected to the main body in advance. For safe use, wear it on the wrist before use. To reconnect the strap after removing it, connect to the bottom of the main body tightly. (If it is difficult to connect, use a pin, etc.)

Caring your Facial Skin

It is recommended to use the product once a week, 300 shots each time.

1 Left cheek (80 shots) → 2 Right cheek (80 shots) → 3 Left forehead (20 shots) → 4 Right forehead (20 shots) → 5 Left temple (10 shots) → 6 Right temple (10 shots) → 7 Left double chin (10 shots) → 8 Center double chin (20 shots) → 9 Right double chin (10 shots) → 10 Left neck (20 shots) → 11 Right neck (20 shots) → End the guide.

※ It is okay to apply shots only to the areas you want, but make sure that you give at least a week of rest to the areas applied with shots.

※ While using the product, the film area of the irradiation part may swell up. As a temporary phenomenon due to the characteristic of product use, it returns to the original state when placed at a room temperature, so use with confidence.

Make sure that you use the product after applying a moisture gel.

If you don't apply a sufficient amount of moisture gel or the gel dries out, shots may not be fired properly or there can be skin irritations such as tingling. Make sure that you apply a sufficient amount of moisture gel to the applicator and the facial skin when you use the product.


For more details on how to use the product, please refer to the Care Guide provided with the product.

Replacing the Cartridge

Check if the power is turned off before replacing it.


※ The cartridge is a consumable. After all shots are used, purchase and replace to a new cartridge.

Separating the cartridge from the main body

- Hold the product main body with one hand and press the Mount/Unmount button.
- At the same time, push up the lower part of the cartridge with the other hand to separate it from the main body.

Mounting the cartridge on the main body

Place the main body with the top/bottom and insertion areas on both sides properly, and push until you hear a click sound to fit it.

- Check if there are any foreign substances (dust, cosmetics, etc.) inside before connecting the cartridge. If there are foreign substances, use a cotton swab or dry cloth to wipe it off before connecting it.

LG Pra.L

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Care Guide

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Facial Contouring Device

BLQ1

詳細說明可查閱 LG Electronics 官方網站的線上使用手冊。
你可以使用智能手機掃描二維碼，快速查閱使用手冊。
Detailed instructions can be found in the WEB manual from
LG Electronics official website. You can quickly access the
user manual by scanning the QR code with a smart device.

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
Caring your Facial Skin


It is recommended to use the product once a week, 300 shots each time.
① Left cheek (80 shots) → ② Right cheek (80 shots) → ③ Left forehead (20 shots) → ④ Right forehead (20 shots) → ⑤ Left temple (10 shots) → ⑥ Right temple (10 shots) → ⑦ Left double chin (10 shots) → ⑧ Center double chin (20 shots) → ⑨ Right double chin (10 shots) → ⑩ Left neck (20 shots) → ⑪ Right neck (20 shots) → End the guide.
※ It is okay to apply shots only to the areas you want, but make sure that you give at least a week of rest to the areas applied with shots.
※ While using the product, the film area of the irradiation part may swell up. As a temporary phenomenon due to the characteristic of product use, it returns to the original state when placed at a room temperature, so use with confidence.

Make sure that you use the product after applying a moisture gel.
If you don't apply a sufficient amount of moisture gel or the gel dries out, shots may not be fired properly or there can be skin irritations such as tingling. Make sure that you apply a sufficient amount of moisture gel to the applicator and the facial skin when you use the product.
※ For detailed information on product use, refer to the Quick Guide provided together.

How to Use


Wash your face thoroughly before starting the care.

Press and hold the Power button for 2 seconds or longer to turn on the product.

If you press the Power button briefly, the product may not turn on or off. Press and hold the Power button. (When turning on: Press and hold for 2 seconds or longer. When turning off: Press and hold for 3 seconds or longer.)


Press the LEVEL button (intensity) button to select intensity (level 1-5).


- ※ When using for the first time, make sure to start from the intensity level 1.
- If skin improves as use time increases, or skin compliance for the pain increases, you may feel less irritation when using in the same energy continuously. If you do not feel irritation, use it by increasing 1-2 level(s).

Apply moisture gel thoroughly in the irradiation part on the device or desired guide areas.


Wrap the product so that the finger touches the SHOT button.

Press the SHOT button with the finger for use. Each time you press the button, 1 shot is irradiated.

- In case of uncomfortable posture during the care, switch the finger pressing the button or hand holding the product to perform the care.

How to Use the Device Correctly


- During the care, apply force to the device and adhere the irradiation part to skin properly for use.
- Do not irradiate to the part same as the one irradiated previously.
- Place the irradiation part in a proper angle to use the product.
- For parts with thin skin, such as neck, use in level 1.

Parts where Care Must Be Avoided

- Care must be avoided between the eyebrows, around the eyes, around the mouth (1cm from the oral angle), jawbone line, and uvula line.
- Adjust the intensity according to the skin properties and do not fire shots to the same area repeatedly or continuously.
- Adjust the number of shots according to the size of a user's guide area when using the device.
- The manufacturer is not responsible for any side effect occurred by providing care to the part that must be avoided.
- 15 dots are irradiated per shot sequentially as below.

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
Safety cautions

- CAUTION concerning the Power Cord**
The Power Plug is the disconnecting device. In case of an emergency, the Power Plug must remain readily accessible.
- CAUTION:** The DC Adapter must be purchased separately, and a certified DC adapter (5 V, 2.0 A or higher) must be used. If you use an uncertified charger, there is a risk of breakdown and change in charging time. If charging is unsuccessful, it is recommended to try again with an alternative certified DC adapter.
- When charging the product, be cautious not to get the charging part to be stained with foreign substances such as moisture or cosmetics. Otherwise, it may cause damage to the product or a fire.
- For safety precautions, refer to the manual provided separately.

Left cheek (80 shots) → Right cheek (80 shots)


Left cheek Front - Middle - Back (40 shots) / Top - Middle - Bottom (40 shots) → Right cheek Front - Middle - Back (40 shots) / Top - Middle - Bottom (40 shots)

- (Common) Move towards each arrow gradually and shoot 1 shot each time. Make sure to repeat shooting the next 1 shot after shooting 1 shot is completely finished each time. For each arrow, perform care as much as the indicated number.
- After irradiating 1 shot, move the product from bottom to top gradually and shoot 1 shot each time along with the facial contour.
- During the care, maintain the arm angle at about 45 degrees for comfortable use.

Forehead (40 shots)


Left bottom - Top (20 shots) → Right bottom - Top (20 shots)

- After irradiating 1 shot, move the product towards each arrow direction and shoot 1 shot each time along with the forehead curve.
- Place the product vertically as shown in the image, and perform care by moving it outward on your face.
- During the care, maintain the arm angle perpendicular to the face and move the product left and right gradually.

Temple (20 shots)


Left front - Bottom (10 shots) → Right front - Bottom (10 shots)

- After irradiating 1 shot, move the product towards each arrow direction and shoot 1 shot each time along with the temple contour.
- Raise the angle of the arm holding the product and adhere the product to the skin to perform care. If the angle is uncomfortable, change a wrist direction left and right to perform care.
- Raise your arm sufficiently to perform care.
- Since it is a narrow part, the product must be used delicately.

Double Chin (40 shots)


Left back - Inside (10 shots) → Middle (20 shots) → Right back - Inside (10 shots)

- After irradiating 1 shot, move the product towards each arrow direction and shoot 1 shot each time as if pulling it forward.
- During the care, make sure that the irradiation part does not face towards sky by tilting your face to adhere it to the skin.
- Since it is a narrow part, the product must be used delicately.
- In case of uncomfortable posture during the care, switch the finger pressing the button or hand holding the product to perform the care.

Neck (40 shots)

Left back - Inside (20 shots) → Right back - Inside (20 shots)

- After irradiating 1 shot, move the product towards each arrow direction and shoot 1 shot each time as if pulling it upward.
- Since it is a narrow part, the product must be used delicately.
- In case of uncomfortable posture during the care, switch the finger pressing the button or hand holding the product to perform the care.